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Method	Performance characteristic	Implementation verification	(food) item verification
Qualitative	Estimated LOD <sub>50</sub> (eLOD <sub>50</sub> )	✓	✓
Quantitative	Intralaboratory reproducibility (S <sub>IR</sub> )	✓	Not applicable
	Estimated bias (eBias)	Not applicable	✓

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Method	Performance characteristic	Implementation verification	(food) item verification
Qualitative	Estimated LOD <sub>50</sub> (eLOD <sub>50</sub> )	Not applicable	✓
Quantitative	Intralaboratory reproducibility (S <sub>IR</sub> )	Not applicable	Not applicable
	Estimated bias (eBias)	Not applicable	✓

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9. Bartram J. and Balance R. World Health Organization. Water Quality Monitoring – A practical to the Design and Implementation of Freshwater Quality Studies and Monitoring Programmes. 1996. UNEP/WHO : Geneva. Available from : [https://www.who.int/water\\_sanitation\\_health/resourcesquality/wqmchap10.pdf](https://www.who.int/water_sanitation_health/resourcesquality/wqmchap10.pdf)
10. Public Health Ontario. Drinking Water Quality – Indicator Bacteria (Total Coliforms and *E. coli*). 2020. Available from : <http://www.publichealthontario.ca>
11. World Health Organization. Guidelines for drinking-water quality. 4<sup>th</sup> ed. Geneva; 2017. p. 297 -297.
12. Ashbolt, N.J., Grabow, W.O.K., Snozzi, M., 2001. Indicators of microbial water quality. In: Fewtrell, L., Bartram, J. (Eds.), Water Quality: Guidelines, Standards and Health. Risk assessment and management for water-related infectious disease. IWA Publishing, London (Chapter 13), pp. 289–315.
13. Leclerc H., Mossel D.A.A., Edberg S.C., and Struijk C.B. Advances in the Bacteriology of the Coliform Group: Their Suitability as Markers of Microbial Water Safety. Annual Review of Microbiology. [Internet]; 2001. 55 : 201-234. Available from : <https://www.annualreviews.org/doi/10.1146/annurev.micro.55.1.201>
14. New York State Department of Health. Coliform Bacteria in Drinking Water Supplies. 2017. [Internet]. Available from : [https://www.health.ny.gov/environmental/water/drinking/coliform\\_bacteria.htm](https://www.health.ny.gov/environmental/water/drinking/coliform_bacteria.htm)
15. ISO/IEC 17025:2017. General requirements for the competence of testing and calibration laboratories. Geneva: International Organization for Standardization; 2017.

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ตาม Standard Method for the Examination of Water and Wastewater 23rd Edition : 2017

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**9221 B. Standard Total Coliform Fermentation Technique**

**1. Samples**

Collect samples as directed in Section 9060A, using sample containers specified in Section 9030B.19. Follow the QC guidelines for sample bottles described in Section 9020B.5*d*. Ensure that samples meet laboratory-acceptance criteria upon receipt.

**2. Quality Control**

All phases of the fermentation technique (9221B–G) require adherence to the quality assurance/quality control (QA/QC) guidelines presented in Section 9020, including, but not limited to, analytical QC (Section 9020B.9), instrumentation/equipment (Sections 9020B.4 and 9030B), and supplies (Section 9020B.5). Refer to Table 9020:I for key QC procedures. Also, note the sections pertaining to appropriate storage and preparation of dehydrated culture media and water quality (Sections 9050 and 9020B.5*f*).

Use commercial dehydrated media when possible, and ensure that their formulations match those specified here because commercial formulations may vary. Prepared fermentation media can be stored in tightly capped tubes or bottles for up to 3 months in the dark, if temperatures are between 1 and 30°C and evaporation is less than 10% of the original volume. If the tubes were refrigerated after sterilization, they should be incubated overnight at room temperature (20°C) before use and those showing growth or bubbles should be discarded to avoid false-positive results. To demonstrate acceptable medium performance, positive and negative culture controls should be tested before first use and as otherwise specified (see Table 9020:VI). Sterility, volume per tube, and pH should also be verified and recorded. To demonstrate comparability between batches of media, perform a use test [Section 9020B.5*f*2)].

If a laboratory is switching to the multiple-tube fermentation technique, analysts ideally should first conduct parallel tests with the previous method to demonstrate applicability and comparability. The results of many coliform performance studies are available in the literature, and the rates of false-positive and -negative results can differ among various media. Users should carefully select the medium and procedure that best fits their needs.

### 3. Presumptive Phase

Use lauryl tryptose broth in this phase of the multiple-tube test, following the QC guidelines cited in 9221B.2.

**a. Reagents and culture medium: Lauryl tryptose broth:**

Tryptose. . . . .	20.0 g
Lactose . . . . .	5.0 g
Dipotassium hydrogen phosphate (K <sub>2</sub> HPO <sub>4</sub> ). . . . .	2.75 g
Potassium dihydrogen phosphate (KH <sub>2</sub> PO <sub>4</sub> ) . . . . .	2.75 g
Sodium chloride (NaCl). . . . .	5.0 g
Sodium lauryl sulfate. . . . .	0.1 g
Reagent-grade water . . . . .	1 L

Add dehydrated ingredients to water, mix thoroughly, and heat to dissolve. Before sterilization, dispense enough medium into fermentation tubes containing inverted vials (also known as Durham tubes) to cover the inverted vial at least one-half to two-thirds after sterilization. Alternatively, omit the inverted vial and add 0.01 g/L bromocresol purple to lauryl tryptose broth (to determine acid production, an indicator of a positive result in this part of the coliform test). Close tubes with metal or heat-resistant plastic caps.

Prepare in accordance with Table 9221:I, making lauryl tryptose broth concentrated enough that adding 100-, 20-, or 10-mL portions of sample to the medium will not reduce ingredient concentrations below those of the standard medium. Autoclave medium at 121°C for 12 to 15 min. Ensure that inverted vials, if used, are free of air bubbles. Medium pH should be 6.8±0.2 after sterilization.

TABLE 9221:I. PREPARATION OF LAURYL TRYPTOSE BROTH

Inoculum <i>mL</i>	Amount of Medium in Tube <i>mL</i>	Volume of Medium Inoculum <i>mL</i>	Dehydrated Lauryl Tryptose Broth Required <i>g/L</i>
1	10 or more	11 or more	35.6
10	10	20	71.2
10	20	30	53.4
20	10	30	106.8
100	50	150	106.8
100	35	135	137.1
100	20	120	213.6

**b. Procedure:**

1) Arrange fermentation tubes in rows of five or ten tubes each in a test tube rack. The number of rows and the sample volumes selected depend on the quality and character of the water to be examined. For potable water, 100 mL must be tested. Use five 20-mL portions, ten 10-mL portions, or one 100-mL portion (a single bottle). For nonpotable water, use five tubes per dilution (of 10, 1, 0.1 mL, etc.).

When making dilutions and measuring diluted sample volumes, follow the precautions given in Section 9215B.2. Use Figure 9215:1 as a guide to preparing dilutions. Shake sample and dilutions vigorously 5 s (about 25 times). Inoculate each tube in a set of five with replicate sample volumes in increasing decimal dilutions, if decimal quantities of the sample are used. Mix test portions in the medium by gentle agitation.

2) Promptly incubate inoculated tubes or bottles, any culture controls, and/or sterility blanks at 35.0°C. After 24 h, swirl each tube or bottle gently and examine it for growth, gas, and/or acidic reaction (shades of yellow color) and, if no gas or acidic reaction is evident, re-incubate and re-examine at the end of 48 h. Record presence or absence of growth, gas, and/or acid production. If the inner vial is omitted, growth with acidity (yellow color) signifies a presumptive-positive reaction.

**c. Interpretation:** Detection of an acidic reaction (yellow color) and/or gas in the tubes or bottles within 48 h constitutes a presumptive-positive reaction. Submit tubes or bottles with a presumptive-positive reaction to the confirmed phase (9221B.4).

The absence of acidic reaction and/or gas formation at the end of 48 h of incubation constitutes a negative test. Submit drinking water samples demonstrating growth without a positive gas or acidic reaction to the confirmed phase (9221B.4).

**4. Confirmed Phase**

**a. Culture medium:** Use BGLB broth fermentation tubes for the confirmed phase, following QC guidelines cited in 9221B.2.

*Brilliant green lactose bile broth:*

Peptone . . . . .	10.0g
Lactose. . . . .	10.0 g
Oxgall . . . . .	20.0 g
Brilliant green. . . . .	0.0133 g
Reagent-grade water. . . . .	1 L

Add dehydrated ingredients to water, mix thoroughly, and heat to dissolve. Before sterilization, dispense medium into fermentation tubes with an inverted vial, ensuring sufficient volume of medium to cover the inverted vial at least one-half to two-thirds after sterilization. Close tubes with metal or heat-resistant plastic caps. Autoclave medium at 121°C for 12 to 15 min. Ensure that inverted vials are free of air bubbles. Medium pH should be  $7.2\pm 0.2$  after sterilization.

**b. Procedure:** Promptly submit all presumptive tubes or bottles showing growth, any amount of gas, or acidic reaction within  $24\pm 2$  h of incubation to the confirmed phase. If additional presumptive tubes or bottles show active fermentation or acidic reaction at the end of a  $48\pm 3$  h incubation period, promptly submit these to the confirmed phase. To confirm presumptive coliform colonies growing on a solid medium using fermentation media, see Section 9222B.4g.

Gently shake or rotate presumptive tubes or bottles showing gas or acidic growth to resuspend the organisms. With a sterile loop 3.0 to 3.5 mm in diameter, transfer one or more loopfuls of culture to a fermentation tube containing BGLB broth. Alternatively, insert a sterile wooden applicator at least 2.5 cm into the culture, promptly remove, and plunge applicator to the bottom of fermentation tube containing BGLB broth. Remove and discard applicator. Repeat for all other presumptive-positive tubes. Analysts may simultaneously inoculate BGLB broth for total coliforms and EC broth for thermotolerant (fecal) coliforms (see 9221E) or EC-MUG broth for *Escherichia coli* (see 9221F). However, if using the same loop or wooden applicator stick to inoculate a culture into more than one medium, inoculate the most inhibitory medium (BGLB broth) last.

Promptly incubate the inoculated BGLB broth tubes at  $35\pm 0.5^\circ\text{C}$ . Any amount of gas formed in the inverted vial of the BGLB broth fermentation tube at any time within  $48\pm 3$  h constitutes a positive confirmed phase. To estimate the coliform density, calculate the MPN value from the number of positive BGLB tubes as described in 9221C.

**c. Alternative procedure:** Use this alternative only for polluted water or wastewater known to produce positive results consistently.

If all presumptive tubes are positive in two or more consecutive dilutions within 24 h, then only submit to the confirmed phase the highest-dilution tubes (smallest sample inoculum) in which all tubes are positive, along with any positive tubes in still higher

dilutions. Submit to the confirmed phase all tubes in which gas or acidic growth is produced in 24 to 48 h.

## 5. Completed Phase

The completed test as described here is not required for drinking-water compliance sample analyses. For nonpotable water samples collected under the Clean Water Act, the requirement that 10% of all total-coliform-positive tubes be subjected to the completed test on a seasonal basis no longer exists. The completed test is included here as a QC recommendation and for use when testing results are uncertain. As additional testing for thermotolerant (fecal) coliforms and/or *E. coli* is required of positive coliform tests, further testing using EC and/or EC-MUG broths is considered a completed test. For QC purposes, if no positive drinking-water samples are received within a quarter, then analyze at least one positive source-water sample to confirm that media respond appropriately.

To verify the presence of coliform bacteria and to provide QC data for nonpotable water-sample analysis, use the completed test on at least one positive sample per quarter. If no positive sample occurs within a quarter, perform a QC check using a known positive sample. Analysts may simultaneously inoculate presumptive-positive media into both BGLB broth for confirmation of total coliforms and EC broth for thermotolerant (fecal) coliforms (9221E) or EC MUG broth for *Escherichia coli* (9221F) as long as BGLB broth is inoculated last. Positive results from incubation in EC and/or EC-MUG broths at elevated temperature ( $44.5\pm 0.2^{\circ}\text{C}$ ) can be considered a completed test. Parallel positive BGLB broth cultures with negative EC or EC-MUG broth cultures indicate the presence of nonfecal coliforms. Parallel positive EC or EC-MUG tubes and negative BGLB broth cultures indicate the presence of thermotolerant (fecal) coliforms or *E. coli*, respectively. Alternatively, the completed test for positive total coliforms may be performed as follows.

**a. Culture media and reagents:** Follow the QC guidelines cited in 9221B.2.

1) *LES Endo agar*—See Section 9222B.2a. Use 100-15-mm Petri plates.

2) *MacConkey agar*:

Peptone . . . . .	17 g
Proteose peptone. . . . .	3 g
Lactose . . . . .	10 g
Bile salts . . . . .	1.5 g
Sodium chloride (NaCl) . . . . .	5 g

Agar . . . . .	13.5 g
Neutral red. . . . .	0.03 g
Crystal violet . . . . .	0.001 g
Reagent-grade water . . . . .	1 L

Add ingredients to water, mix thoroughly, and heat to boiling to dissolve. Sterilize by autoclaving for 15 min at 121°C. Temper agar after sterilization and pour into Petri plates (100×15 mm). Medium pH should be 7.1±0.2 after sterilization.

3) *Nutrient agar*:

Peptone . . . . .	5.0 g
Beef extract . . . . .	3.0 g
Agar . . . . .	15.0 g
Reagent-grade water . . . . .	1 L

Add ingredients to water, mix thoroughly, and heat to dissolve. Before sterilization, dispense in screw-capped tubes. Autoclave at 121°C for 15 min. Medium pH should be 6.8 ±0.2 after sterilization. After sterilization, immediately place tubes in an inclined position so the agar will solidify with a sloped surface. Tighten screw caps after cooling and store in a protected, cool storage area.

4) *Gram-stain reagents*—Reagents are commercially available as prepared solutions.

a) *Ammonium oxalate-crystal violet (Hucker's)*—Dissolve 2 g crystal violet (90% dye content) in 20 mL 95% ethyl alcohol. **CAUTION: Flammable.** Dissolve 0.8 g  $(\text{NH}_4)_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$  in 80 mL reagent-grade water. Mix the two solutions and age for 24 h before use. Filter through paper into a staining bottle.

b) *Lugol's solution, Gram's modification*—Grind 1 g iodine crystals and 2 g KI in a mortar. Add reagent-grade water, a few milliliters at a time, and grind thoroughly after each addition until solution is complete. Rinse solution into an amber glass bottle with the remaining water, using a total of 300 mL.

c) *Counterstain*—Dissolve 2.5 g safranin dye in 100 mL 95% ethyl alcohol. Add 10 mL to 100 mL reagent-grade water. **CAUTION: Flammable.**

d) *Acetone alcohol*—Mix equal volumes of ethyl alcohol (95%) with acetone. **CAUTION: Flammable.**

**b. Procedure:**

1) Using aseptic technique, streak one LES Endo agar (Section 9222B.2a) or MacConkey agar plate from each presumptive positive tube of BGLB broth as soon as possible after gas is observed. Streak plates in a manner to ensure the presence of some discrete colonies separated by at least 0.5 cm. To obtain a high proportion of successful isolations if coliform organisms are present, use the following approach:

- a) Use a sterile 3-mm-diam loop or an inoculating needle slightly curved at the tip;
- b) tap and incline the fermentation tube to avoid picking up any membrane or scum on the needle;
- c) insert the end of the loop or needle into the liquid in the tube to a depth of approximately 0.5 cm; and
- d) streak a plate for isolation with the curved section of the needle in contact with the agar to avoid a scratched or torn surface. Flame the loop between the second and third quadrants to improve colony isolation. Incubate plates, inverted, at  $35 \pm 0.5^\circ\text{C}$  for  $24 \pm 2$  h.

2) The colonies developing on LES Endo agar are defined as *typical* (pink to dark red with a green metallic surface sheen) or *atypical* (pink, red, white, or colorless colonies without sheen) after 24 h incubation. Typical lactose-fermenting colonies developing on MacConkey agar are red and may be surrounded by an opaque zone of precipitated bile. From each plate, pick one or more typical, well-isolated coliform colonies or, if no typical colonies are present, pick two or more colonies considered most likely to be coliforms. Transfer growth from each isolate to a single-strength lauryl tryptose broth fermentation tube and onto a nutrient agar slant.

If needed, use a colony-magnifying device to provide optimum magnification when colonies are picked from the LES Endo or MacConkey agar plates. When transferring colonies, choose well-isolated ones and barely touch the colony surface with a flame-sterilized, air-cooled transfer needle to minimize the danger of transferring a mixed culture. Incubate secondary broth tubes (lauryl tryptose broth with inverted fermentation vials) at  $35 \pm 0.5^\circ\text{C}$  for  $24 \pm 2$  h; if gas is not produced within  $24 \pm 2$  h, reincubate and examine again at  $48 \pm 3$  h. Microscopically examine Gram-stained preparations from those 24-h nutrient agar slant cultures corresponding to the secondary tubes that show gas.

3) Gram-stain technique—The Gram stain may be omitted from the completed test for potable-water samples only because Gram-positive bacteria and spore-forming organisms in drinking water rarely survive this selective screening procedure.

Various modifications of the Gram stain technique exist. Use Hucker's modification (as follows) for staining smears of pure cultures; include a Gram-positive and a Gram-negative culture as controls.

On one slide, prepare separate light emulsions of the test bacterial growth and positive and negative control cultures using drops of distilled water on the slide. Air-dry, fix by passing slide through a flame, and stain for 1 min with ammonium oxalatecrystal violet solution. Rinse slide in tap water and drain off excess; apply Lugol's solution for 1 min.

Rinse stained slide in tap water. Decolorize for approximately 15 to 30 s with acetone alcohol by holding the slide between the fingers and letting acetone alcohol flow across the stained smear until the solvent flows colorlessly from the slide. Do not overdecolorize. Counterstain with safranin for 15 s, rinse with tap water, blot dry with absorbent paper or air dry, and examine microscopically. Gram-positive organisms are blue; Gramnegative organisms are red. Results are acceptable only when controls have given proper reactions.

**c. Interpretation:** Formation of gas in the secondary tube of lauryl tryptose broth within  $48 \pm 3$  h and demonstration of Gram-negative, nonspore-forming, rod-shaped bacteria from the agar culture constitute a positive result for the completed test, demonstrating that a member of the coliform group is present.

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**9221 E. Thermotolerant (Fecal) Coliform Procedure**

Traditionally called *fecal coliforms*, thermotolerant coliforms (those that ferment lactose to produce gas at 44.5°C) have been documented in organically rich waters or tropical climates in the absence of recent fecal contamination. So when looking for evidence of fecal contamination, testing for *E. coli*—a more specific indicator—is recommended. Nevertheless, regulations may require that thermotolerant (fecal) coliforms be identified and enumerated.

To test for thermotolerant coliforms, use one of the multiple tube procedures described here or the membrane-filter methods described in Sections 9222D and E. In the multiple-tube fermentation technique, thermotolerant coliforms are identified by their ability to ferment lactose to produce gas at 44.5±0.2°C within 24±2 h.

**1. Thermotolerant Coliform Test (EC Medium)**

The thermotolerant coliform test using EC medium is applicable to investigations of drinking water, stream pollution, unfiltered raw water sources, wastewater treatment systems, bathing waters, seawaters, and general water-quality monitoring. Do not use EC medium to directly isolate thermotolerant coliforms from water; prior enrichment in a presumptive medium is required for optimum recovery of thermotolerant coliforms. (To test presumptive coliform colonies growing on solid media, refer to Section 9222G.3c)

**a. EC medium:** Prepare EC medium following QC guidelines cited in 9221B.2.

Tryptose or trypticase . . . . .	20.0 g
Lactose . . . . .	5.0 g
Bile salts mixture or bile salts No. 3 . . . . .	1.5 g
Dipotassium hydrogen phosphate (K <sub>2</sub> HPO <sub>4</sub> ). . . . .	4.0 g
Potassium dihydrogen phosphate (KH <sub>2</sub> PO <sub>4</sub> ) . . . . .	1.5 g
Sodium chloride (NaCl). . . . .	5.0 g
Reagent-grade water. . . . .	1 L

Add dehydrated ingredients to water, mix thoroughly, and heat to dissolve. Before sterilization, dispense sufficient medium in fermentation tubes with an inverted vial to cover

the inverted vial at least one-half to two-thirds after sterilization. Close tubes with metal or heat-resistant plastic caps. Autoclave medium at 121°C for 12 to 15 min. Ensure that inverted vials are free of air bubbles. Medium pH should be  $6.9 \pm 0.2$  after sterilization.

**b. Procedure:**

1) After incubation, gently shake or rotate fermentation tubes or bottles showing gas, growth, or acidity to resuspend the organisms. Promptly use a sterile 3- to 3.5-mm-diam loop to transfer one or more loopfuls of culture from bottles or tubes showing growth with acid and/or gas production to a fermentation tube containing EC broth. Alternatively, insert a sterile wooden applicator at least 2.5 cm into the culture, promptly remove, and plunge applicator to the bottom of a fermentation tube containing EC broth. Remove and discard applicator. Repeat for all other presumptive-positive tubes and incubate at  $44.5 \pm 0.2^\circ\text{C}$ .

Simultaneous inoculation into EC broth and/or EC-MUG broth along with BGLB broth is acceptable, if the most inhibitory medium (BGLB broth) is inoculated last.

2) Place all EC tubes into a circulating water bath (preferably with a gabled cover) within 30 min after inoculation. Incubate inoculated EC broth tubes at  $44.5 \pm 0.2^\circ\text{C}$  for  $24 \pm 2$  h. Maintain a sufficient water depth in the water bath incubator to immerse tubes to the upper level of the medium.

**c. Interpretation:** Gas production with growth in an EC broth culture within  $24 \pm 2$  h or less is considered a positive thermotolerant (fecal) coliform reaction. Failure to produce gas (with little or no growth) constitutes a negative reaction. If multiple tubes are used, calculate the MPN of thermotolerant coliforms from the number of positive EC broth tubes, as described in 9221C. When using only one tube for subculturing from a single presumptive bottle, report as the presence or absence of thermotolerant coliforms. If heavy growth occurs with no gas production, subject the culture to a thermotolerant coliform or *E. coli* test using a different medium.

## 2. Thermotolerant (Fecal) Coliform Direct Test (A-1 Medium)

**a. A-1 medium:** This medium may be used to directly isolate thermotolerant coliforms from unfiltered source water, treated wastewater, and seawater, but not drinking water. Follow guidelines in 9221B.1 for sample collection. Unlike EC medium, A-1 medium does not require prior enrichment in a presumptive medium for optimum recovery of thermotolerant coliforms. Use QC guidelines cited in 9221B.2.

Lactose . . . . .	5.0 g
Tryptone . . . . .	20.0 g
Sodium chloride (NaCl). . . . .	5.0 g
Salicin. . . . .	0.5 g
Polyethylene glycol p-isooctylphenyl ether* . . . . .	1.0 mL
Reagent-grade water . . . . .	1 L

Heat to dissolve solid ingredients, add polyethylene glycol p-isooctylphenyl ether, and adjust to pH  $6.9 \pm 0.1$ . For 10-mL samples, prepare double-strength medium so the final concentration of ingredients after sample addition is correct. Before sterilization, dispense sufficient medium in fermentation tubes with an inverted vial to cover the inverted vial at least one-half to two-thirds after sterilization. Close with metal or heat-resistant plastic caps. Sterilize by autoclaving at  $121^{\circ}\text{C}$  for 10 min. Ensure that inverted vials are free of air bubbles. Store in the dark at room temperature for not longer than 7 d. Ignore precipitate formed during storage.

**b. Procedure:** Inoculate tubes of A-1 broth as directed in 9221B.3b. Incubate for 3 h at  $35 \pm 0.5^{\circ}\text{C}$ . Transfer tubes to a water bath at  $44.5 \pm 0.2^{\circ}\text{C}$  and incubate for another  $21 \pm 2$  h.

Interpretation: Gas production in any A-1 broth culture within 24 h or less is a positive reaction [i.e., thermotolerant (fecal) coliforms are present]. Calculate the MPN of thermotolerant (fecal) coliforms from the number of positive A-1 broth tubes, as described in 9221C.

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**9221 F. Escherichia coli Procedure Using Fluorogenic Substrate**

*Escherichia coli* is a member of the indigenous fecal flora of warm-blooded animals. The presence of *E. coli* in water is considered a specific indicator of fecal contamination and the possible presence of enteric pathogens. Tests for *E. coli* are applicable to the analysis of drinking, surface, ground, and waste water. Testing for *E. coli* can be performed using the multiple tube procedure described here, by the membrane filter method described in Section 9222G, or by the chromogenic enzyme substrate tests described in Section 9223. Other *E. coli* procedures are presented in 9221G.

For the *E. coli* test using EC-MUG medium, *E. coli* is defined as the species of coliform bacteria that possesses the enzyme  $\beta$ -glucuronidase, which can cleave the fluorogenic substrate 4-methylumbelliferyl-  $\beta$ -D-glucuronide (MUG), thus releasing the fluorogen within  $24 \pm 2$  h or less when grown in EC-MUG medium at  $44.5 \pm 0.2^\circ\text{C}$ .

**1. Escherichia coli Test (EC-MUG Medium)**

The use of EC-MUG medium to detect *E. coli* is applicable to investigations of drinking water, stream pollution, unfiltered raw water sources, wastewater treatment systems, bathing waters, seawaters, and general water-quality monitoring. Do not use EC-MUG for the direct isolation of *E. coli*; prior enrichment in a presumptive medium is required for optimum recovery. (To test presumptive coliform colonies growing on solid media, refer to Section 9222G.2.)

Use EC-MUG medium to test for *E. coli* in a total coliform positive culture, following QC guidelines cited in 9221B.2.

**a. EC-MUG medium:** Prepare EC-MUG medium following QC guidelines cited in 9221B.2.

Tryptose or trypticase . . . . .	20.0 g
Lactose . . . . .	5.0 g
Bile salts mixture or bile salts No. 3. . . . .	1.5 g
Dipotassium hydrogen phosphate (K <sub>2</sub> HPO <sub>4</sub> ). . . . .	4.0 g
Potassium dihydrogen phosphate (KH <sub>2</sub> PO <sub>4</sub> ) . . . . .	1.5 g

Sodium chloride (NaCl) . . . . .	5.0 g
4-Methylumbelliferyl- -D-glucuronide (MUG) . . . . .	0.05 g
Reagent-grade water . . . . .	1 L

Add dehydrated ingredients to water, mix thoroughly, and heat to dissolve. Before sterilization, dispense in tubes that do not fluoresce under long-wavelength (365–366 nm) ultraviolet (UV) light. An inverted tube is not necessary. Close tubes with metal or heat-resistant plastic caps. Medium pH should be  $6.9 \pm 0.2$  after sterilization for 15 min at  $121^\circ\text{C}$ .

**b. Procedure:**

1) Gently shake or rotate fermentation tubes or bottles showing growth, gas, or acidity to resuspend the organisms. Using a sterile 3- or 3.5-mm-diam loop, transfer one or more loopfuls of growth from the fermentation tube or bottle to EC-MUG broth. Alternatively, insert a sterile wooden applicator stick at least 2.5 cm into the culture, promptly remove, and plunge applicator to the bottom of a fermentation tube containing EC-MUG broth.

2) Place all EC-MUG tubes in water bath within 30 min after inoculation. Incubate inoculated EC-MUG tubes and negative controls for  $24 \pm 2$  h in a circulating water bath (preferably with a gable cover) maintained at  $44.5 \pm 0.2^\circ\text{C}$ . Maintain a sufficient water depth in the water-bath incubator to immerse tubes to the upper level of medium.

**c. Interpretation:** Examine all tubes exhibiting growth for fluorescence using a 6W, 365–366 nm long-wavelength UV lamp. The presence of bright blue fluorescence is considered a positive result for *E. coli*. Growth in the absence of bright blue fluorescence is considered a negative result. To help interpret results and avoid misidentifying weak autofluorescence of the medium or glass tubes as a positive response, include in the assay a positive control [a known *E. coli* (MUG-positive) culture], a negative control [a thermotolerant *Klebsiella pneumoniae* (MUG-negative) culture], and an uninoculated medium control. The distance between the UV lamp and the tubes should be such that the *E. coli* positive control shows distinct fluorescence while the MUG-negative and uninoculated controls do not. If using multiple tubes, calculate the MPN for *E. coli* from the number of positive EC-MUG broth tubes, as described in 9221C. When using only one tube, or subculturing from a single presumptive bottle or colony, report as the presence or absence of *E. coli*.

## 2. Simultaneous Determination of Thermotolerant Coliforms and *E. coli*

The presence of thermotolerant coliforms and *E. coli* can be determined simultaneously by including an inverted vial (Durham tube) in tubes of EC-MUG broth. Prepare EC-MUG broth according to 9221F.1.

**a. Setup:** Before sterilization dispense, in fermentation tubes with an inverted vial, sufficient medium to cover the inverted vial at least one-half to two-thirds after sterilization. Close with metal or heat-resistant caps. Medium pH should be  $6.9 \pm 0.2$  after sterilization for 15 min at  $121^\circ\text{C}$ .

**b. Procedure:**

1) Gently shake or rotate fermentation tubes or bottles showing growth, gas, or acidity to resuspend the organisms. Using a sterile 3- or 3.5-mm-diam loop, transfer one or more loopfuls of growth from the fermentation tube or bottle to EC-MUG broth. Alternatively, insert a sterile wooden applicator stick at least 2.5 cm into the culture, promptly remove, and plunge applicator to the bottom of a fermentation tube containing EC-MUG broth.

2) Place all EC-MUG tubes in water bath within 30 min after inoculation. Incubate inoculated EC-MUG tubes, along with positive and negative controls, for  $24 \pm 2$  h in a circulating water bath (preferably with a gable cover) maintained at  $44.5 \pm 0.2^\circ\text{C}$ . Maintain a sufficient water depth in the water-bath incubator to immerse tubes to the upper level of medium.

**c. Interpretation:** Examine all tubes exhibiting growth and/or gas for fluorescence using a 6W, 365–366 nm long-wavelength UV lamp. Growth with gas production is considered a positive result for thermotolerant coliforms. The presence of bright blue fluorescence is considered a positive result for *E. coli*. Tubes with growth and/or gas *and* fluorescence are considered positive for both thermotolerant coliforms and *E. coli*. Tubes with growth and/or gas but without bright blue fluorescence are considered positive for thermotolerant coliforms and negative for *E. coli*.

Due to indigenous autofluorescence of media or glass tubes/ inserts, use caution in interpreting results. To help interpret results, include in each assay a positive control [a known *E. coli* (MUG-positive) culture], a negative control [a thermotolerant *Klebsiella*

*pneumoniae* (MUG-negative) culture], and an uninoculated medium control. The distance between the UV lamp and the tubes should be such that the *E. coli* positive control shows distinct fluorescence while the MUG-negative and uninoculated controls do not. If multiple tubes are used, calculate the MPN for *E. coli* and thermotolerant coliforms from the number of positive EC-MUG broth tubes, as described in 9221C. When using only one tube, or subculturing from a single presumptive bottle or colony, report the presence or absence of *E. coli* and thermotolerant coliforms.